

Mob3B rabbit pAb

Cat No.:ES2811

For research use only

Overview

Product Name	Mob3B rabbit pAb
Host species	Rabbit
Applications	WB;IHC;IF;ELISA
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	Western Blot: 1/500 - 1/2000.
	Immunohistochemistry: 1/100 - 1/300. ELISA:
	1/5000. Not yet tested in other applications.
Immunogen	The antiserum was produced against synthesized
0	peptide derived from human MOBKL2B. AA
	range:71-120
Specificity	Mob3B Polyclonal Antibody detects endogenous
	levels of Mob3B protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	MOB kinase activator 3B
Gene Name	MOB3B
Cellular localization	
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	23kD
Human Gene ID	79817
Human Swiss-Prot Number	Q86TA1
Alternative Names	MOB3B; MOBKL2B; MOB kinase activator 3B; Mob1
	homolog 2b; Mps one binder kinase activator-like 2B
Background	The protein encoded by this gene shares similarity
	with the yeast Mob1 protein. Yeast Mob1 binds
	Mps1p, a protein kinase essential for spindle pole
	body duplication and mitotic checkpoint regulation.
	This gene is located on the opposite strand as the



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



interferon kappa precursor (IFNK) gene. [provided by RefSeq, Jul 2008],

Western Blot analysis of various cells using Mob3B Polyclonal Antibody

Immunohistochemical analysis of paraffin-embedded Human colon cancer. Antibody was diluted at 1:100(4° overnight). High-pressure and temperature Tris-EDTA,pH8.0 was used for antigen retrieval. Negetive contrl (right) obtaned from antibody was pre-absorb



COLO205

(kD)

117-85-

48-

34-

26-

19-

-- 117 -- 85 -- 48 -- 34 -- 26 -- 19 (kD) Western blot analysis of lysates from COLO cells, using MOBKL2B Antibody. The lane on the right is blocked with the synthesized peptide.



www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C