

## MRP-S25 rabbit pAb

Cat No.: ES2849

For research use only

## Overview

Product Name MRP-S25 rabbit pAb

Host species Rabbit

Applications WB;ELISA;IHC Species Cross-Reactivity Human;Mouse;Rat

**Recommended dilutions** WB 1:500-2000;IHC-p 1:50-300; ELISA 2000-20000 **Immunogen** The antiserum was produced against synthesized

peptide derived from human MRPS25. AA

range:124-173

**Specificity** MRP-S25 Polyclonal Antibody detects endogenous

levels of MRP-S25 protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

**Storage** Store at -20°C. Avoid repeated freeze-thaw cycles.

**Protein Name** 28S ribosomal protein S25 mitochondrial

Gene Name MRPS25

Cellular localization Mitochondrion .

**Purification** The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 20kD
Human Gene ID 64432
Human Swiss-Prot Number P82663

Alternative Names MRPS25; RPMS25; 28S ribosomal protein S25;

mitochondrial; MRP-S25; S25mt

**Background** Mammalian mitochondrial ribosomal proteins are

encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition

compared to prokaryotic ribosomes, where this ratio




+86-27-59760950 ELKbio@ELKbiotech.com

www.elkbiotech.com




is reversed. Another difference between mammalian mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 28S subunit protein. A pseudogene corresponding to this gene is found on chromosome 4. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2016],

Western Blot analysis of various cells using MRP-S25 Polyclonal Antibody diluted at 1:2000




Western blot analysis of lysates from 293 cells, using MRPS25 Antibody. The lane on the right is blocked with the synthesized peptide.




+86-27-59760950







Western blot analysis of the lysates from HepG2 cells using MRPS25 antibody.



Immunohistochemical analysis of paraffin-embedded human Squamous cell carcinoma of lung. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).

