

## MRP-S9 rabbit pAb

Cat No.: ES2854

For research use only

## Overview

Product Name MRP-S9 rabbit pAb

Host species Rabbit

Applications WB;IHC;IF;ELISA Species Cross-Reactivity Human;Mouse;Rat

**Recommended dilutions** Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.

Immunogen The antiserum was produced against synthesized

peptide derived from human MRPS9. AA

range:91-140

**Specificity** MRP-S9 Polyclonal Antibody detects endogenous

levels of MRP-S9 protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

**Store at -20°C.** Avoid repeated freeze-thaw cycles.

**Protein Name** 28S ribosomal protein S9 mitochondrial

Gene Name MRPS9

Cellular localization Mitochondrion.

**Purification** The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 40kD
Human Gene ID 64965
Human Swiss-Prot Number P82933

Alternative Names MRPS9; RPMS9; 28S ribosomal protein S9;

mitochondrial; MRP-S9; S9mt

**Background** Mammalian mitochondrial ribosomal proteins are

encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an

www.elkbiotech.com




+86-27-59760950 ELKbio@ELKbiotech.com




estimated 75% protein to rRNA composition compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 28S subunit protein. [provided by RefSeq, Jul 2008],

Western Blot analysis of various cells using MRP-S9 Polyclonal Antibody




Western blot analysis of lysates from A549 cells, using MRPS9 Antibody. The lane on the right is blocked with the synthesized peptide.



+86-27-59760950







Immunohistochemical analysis of paraffin-embedded human liver cancer. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).



23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei, P.R.C

+86-27-59760950