

Mucin 13 rabbit pAb

Cat No.: ES2868

For research use only

Overview

Product Name Mucin 13 rabbit pAb

Host species Rabbit

Applications WB;IHC;IF;ELISA Species Cross-Reactivity Human;Rat;Mouse;

Recommended dilutions Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications. The antiserum was produced against synthesized

Immunogen The antiserum was produced against synthesized

peptide derived from human MUC13. AA

range:421-470

Specificity Mucin 13 Polyclonal Antibody detects endogenous

levels of Mucin 13 protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Mucin-13
Gene Name MUC13

Cell ular localization Cell membrane ; Single-pass type I membrane

protein . Apical cell membrane . Secreted . Also

exists as a soluble form.

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 50kD
Human Gene ID 56667
Human Swiss-Prot Number Q9H3R2

Alternative Names MUC13; DRCC1; RECC; Mucin-13; MUC-13;

Down-regulated in colon cancer 1

Background Epithelial mucins, such as MUC13, are a family of

secreted and cell surface glycoproteins expressed by



+86-27-59760950 ELKbio@ELKbiotech.com www.elkbiotech.co



120

90


50

al., 2001 [PubMed 11278439]).[supplied by OMIM,
Jul 2008],

Western Blot analysis of various cells using Mucin 13


Western Blot analysis of various cells using Mucin 13 Polyclonal Antibody diluted at 1:1000


ductal and glandular epithelial tissues (Williams et



Mucin 13


Western Blot analysis of COLO205 cells using Mucin 13 Polyclonal Antibody diluted at 1:1000











Immunohistochemistry analysis of paraffin-embedded human colon carcinoma tissue, using MUC13 Antibody. The picture on the right is blocked with the synthesized peptide.



+86-27-59760950