

Pax-9 rabbit pAb

Cat No.:ES3155

For research use only

Overview

Product Name Pax-9 rabbit pAb

Host species Rabbit
Applications WB;ELISA

Species Cross-Reactivity Human; Mouse; Rat

Recommended dilutions Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not

yet tested in other applications.

Immunogen The antiserum was produced against synthesized

peptide derived from human Pax-9. AA

range:158-207

Specificity Pax-9 Polyclonal Antibody detects endogenous

levels of Pax-9 protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Paired box protein Pax-9

Gene Name PAX9 **Cellular localization** Nucleus.

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 35kD
Human Gene ID 5083
Human Swiss-Prot Number P55771

Alternative Names PAX9; Paired box protein Pax-9

Background This gene is a member of the paired box (PAX) family

of transcription factors. Members of this gene family

typically contain a paired box domain, an octapeptide, and a paired-type homeodomain. These genes play critical roles during fetal

development and cancer growth. Mice lacking this gene exhibit impaired development of organs,




+86-27-59760950 ELKbio@ELKbiotech.com


www.elkbiotech.com



musculature and the skeleton, including absent and abnormally developed teeth, and neonatal lethality. Mutations in the human gene are associated with selective tooth agenesis-3. [provided by RefSeq, Sep 2015],



Western Blot analysis of various cells using Pax-9
Polyclonal Antibody cells nucleus extracted by Minute TM
Cytoplasmic and Nuclear Fractionation kit
(SC-003,Inventbiotech,MN,USA).



+86-27-59760950

Western blot analysis of lysate from rat heart, using Pax-9 antibody.

