

Periphilin 1 rabbit pAb

Cat No.: ES3175

For research use only

Overview

Product Name Periphilin 1 rabbit pAb

Host species Rabbit

Applications WB;IHC;IF;ELISA Species Cross-Reactivity Human;Mouse

Recommended dilutions Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.

Immunogen The antiserum was produced against synthesized

peptide derived from human PPHLN. AA

range:71-120

Specificity Periphilin 1 Polyclonal Antibody detects endogenous

levels of Periphilin 1 protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Periphilin-1
Gene Name PPHLN1

Cellular localization Nucleus . Cytoplasm . Chromosome . In

undifferentiated keratinocytes expressed in speckle-type nuclear granules and at the nuclear membrane, but in the differentiated keratinocytes colocalized with periplakin at the cell periphery and

at cell-cell junctio

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 53kD
Human Gene ID 51535
Human Swiss-Prot Number Q8NEY8

Alternative Names PPHLN1; HSPC206; HSPC232; Periphilin-1; Gastric

cancer antigen Ga50




+86-27-59760950 ELKbio@ELKbiotech.com www.elkbiotech.com




Background

The protein encoded by this gene is one of the several proteins that become sequentially incorporated into the cornified cell envelope during the terminal differentiation of keratinocyte at the outer layers of epidermis. This protein interacts with periplakin, which is known as a precursor of the cornified cell envelope. The cellular localization pattern and insolubility of this protein suggest that it may play a role in epithelial differentiation and contribute to epidermal integrity and barrier formation. Multiple alternatively spliced transcript variants encoding distinct isoforms have been observed. [provided by RefSeq, Jul 2008],

Western Blot analysis of various cells using Periphilin 1 Polyclonal Antibody




Immunohistochemistry analysis of paraffin-embedded human liver carcinoma tissue, using PPHLN Antibody. The picture on the right is blocked with the synthesized peptide.









Western blot analysis of lysates from COLO and 293 cells, using PPHLN Antibody. The lane on the right is blocked with the synthesized peptide.

