

PKA IIβ reg rabbit pAb

Cat No.:ES3217

For research use only

Overview

Product Name	PKA IIβ reg rabbit pAb
Host species	Rabbit
Applications	WB;IHC;IF;ELISA
Species Cross-Reactivity	Human; Mouse; Rat
Recommended dilutions	Western Blot: 1/500 - 1/2000.
	Immunohistochemistry: 1/100 - 1/300. ELISA:
	1/10000. Not yet tested in other applications.
Immunogen	The antiserum was produced against synthesized
	peptide derived from human PKA-R2 beta. AA
	range:79-128
Specificity	PKA IIβ reg Polyclonal Antibody detects endogenous
	levels of PKA IIβ reg protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	cAMP-dependent protein kinase type II-beta
	regulatory subunit
Gene Name	PRKAR2B
Cellular localization	Cytoplasm . Cell membrane . Colocalizes with PJA2 in
	the cytoplasm and at the cell membrane.
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	46kD
Human Gene ID	5577
Human Swiss-Prot Number	P31323
Alternative Names	PRKAR2B; cAMP-dependent protein kinase type
	II-beta regulatory subunit
Background	cAMP is a signaling molecule important for a variety
	of cellular functions. cAMP exerts its effects by
	activating the cAMP-dependent protein kinase,



+86-27-59760950


ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



which transduces the signal through phosphorylation of different target proteins. The inactive kinase holoenzyme is a tetramer composed of two regulatory and two catalytic subunits. cAMP causes the dissociation of the inactive holoenzyme into a dimer of regulatory subunits bound to four cAMP and two free monomeric catalytic subunits. Four different regulatory subunits and three catalytic subunits have been identified in humans. The protein encoded by this gene is one of the regulatory subunits. This subunit can be phosphorylated by the activated catalytic subunit. This subunit has been shown to interact with and suppress the transcriptional activity of the cAMP responsive element binding protein 1 (CREB1) in activ



Western Blot analysis of various cells using PKA II β reg Polyclonal Antibody

Western Blot analysis of 3T3 cells using PKA II β reg Polyclonal Antibody



15----


ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



Immunohistochemistry analysis of paraffin-embedded human colon carcinoma tissue, using PKA-R2 beta Antibody. The picture on the right is blocked with the synthesized peptide.



Western blot analysis of lysates from K562 and HepG2 cells, using PKA-R2 beta Antibody. The lane on the right is blocked with the synthesized peptide.



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C