

Rap 2C rabbit pAb

Cat No.:ES3325

For research use only

Overview

Product Name	Rap 2C rabbit pAb
Host species	Rabbit
Applications	WB;IHC
Species Cross-Reactivity	Human; Mouse
Recommended dilutions	WB 1:500-2000;IHC-p 1:50-300
Immunogen	The antiserum was produced against synthesized
	peptide derived from human RAP2C. AA
	range:97-146
Specificity	Rap 2C Polyclonal Antibody detects endogenous
	levels of Rap 2C protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Ras-related protein Rap-2c
Gene Name	RAP2C
Cellular localization	Cytoplasm . Recycling endosome membrane ;
	Lipid-anchor ; Cytoplasmic side .
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	26kD
Human Gene ID	57826
Human Swiss-Prot Number	Q9Y3L5
Alternative Names	RAP2C; Ras-related protein Rap-2c
Background	The protein encoded by this gene is a member of
	the Ras-related protein subfamily of the Ras GTPase
	superfamily. Members of this family are small
	GTPases that act as molecular switches to regulate
	cellular proliferation, differentiation, and apoptosis.
	This protein has been reported to activate in vitro
	transcriptional activity of the serum response




+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C





element. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2012],

Western Blot analysis of various cells using Rap 2C Polyclonal Antibody diluted at 1:500

Western blot analysis of lysates from Jurkat and MCF-7 cells, using RAP2C Antibody. The lane on the right is blocked with the synthesized peptide.


Western blot analysis of the lysates from K562 cells using RAP2C antibody.



www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C





Immunohistochemical analysis of paraffin-embedded human Breast cancer. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C