

## RASSF4 rabbit pAb

Cat No.: ES3329

For research use only

## Overview

Product Name RASSF4 rabbit pAb

Host speciesRabbitApplicationsWB;ELISASpecies Cross-ReactivityHuman;Rat

Recommended dilutions Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not

yet tested in other applications.

Immunogen The antiserum was produced against synthesized

peptide derived from human RASSF4. AA

range:134-183

**Specificity** RASSF4 Polyclonal Antibody detects endogenous

levels of RASSF4 protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Ras association domain-containing protein 4

Gene Name RASSF4

**Cellular localization** 

**Purification** The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 37kD
Human Gene ID 83937
Human Swiss-Prot Number Q9H2L5

Alternative Names RASSF4; AD037; Ras association domain-containing

protein 4

**Background** The function of this gene has not yet been

determined but may involve a role in tumor


suppression. Alternative splicing of this gene results in several transcript variants; however, most of the variants have not been fully described. [provided by

RefSeq, Jul 2008],




+86-27-59760950 ELKbio@ELKbiotech.com






Western Blot analysis of various cells using RASSF4 Polyclonal Antibody



Western blot analysis of lysates from COLO and HUVEC cells, using RASSF4 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from RAW264.7cells using RASSF4 antibody.

