

Ribosomal Protein L3 rabbit pAb

Cat No.: ES3363

For research use only

Overview

Product Name Ribosomal Protein L3 rabbit pAb

Host species Rabbit

Applications WB;IHC;IF;ELISA Species Cross-Reactivity Human;Mouse;Rat

Recommended dilutions Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.

Immunogen The antiserum was produced against synthesized

peptide derived from human RPL3. AA

range:211-260

Specificity Ribosomal Protein L3 Polyclonal Antibody detects

endogenous levels of Ribosomal Protein L3 protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name 60S ribosomal protein L3

Gene Name RPL3

Cellular localization Nucleus, nucleolus. Cytoplasm.

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 46kD
Human Gene ID 6122
Human Swiss-Prot Number P39023

Alternative Names RPL3; OK/SW-cl.32; 60S ribosomal protein L3; HIV-1

TAR RNA-binding protein B; TARBP-B

Background Ribosomes, the complexes that catalyze protein


synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal




+86-27-59760950 ELKbio@ELKbiotech.com



protein that is a component of the 60S subunit. The protein belongs to the L3P family of ribosomal proteins and it is located in the cytoplasm. The protein can bind to the HIV-1 TAR mRNA, and it has been suggested that the protein contributes to tat-mediated transactivation. This gene is co-transcribed with several small nucleolar RNA genes, which are located in several of this gene's introns. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome. [provided by RefSeq, Jul



Western Blot analysis of various cells using Ribosomal Protein L3 Polyclonal Antibody




+86-27-59760950

Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using RPL3 Antibody. The picture on the right is blocked with the synthesized peptide.







Western blot analysis of lysates from HT-29 cells, using RPL3 Antibody. The lane on the right is blocked with the synthesized peptide.

Western blot analysis of the lysates from HUVECcells using RPL3 antibody.

