

Serpin B12 rabbit pAb

Cat No.:ES3430

For research use only

Overview

Product Name Serpin B12 rabbit pAb

Host species Rabbit
Applications WB;IHC

Species Cross-Reactivity Human; Rat; Mouse;

Recommended dilutions WB 1:500-2000;IHC-p 1:50-300

Immunogen The antiserum was produced against synthesized

peptide derived from human SERPINB12. AA

range:241-290

Specificity Serpin B12 Polyclonal Antibody detects endogenous

levels of Serpin B12 protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein NameSerpin B12Gene NameSERPINB12Cellular localizationCytoplasm .

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 47kD
Human Gene ID 89777
Human Swiss-Prot Number Q96P63

Alternative Names SERPINB12; Serpin B12

Background function:Inhibits trypsin and plasmin, but not

thrombin, coagulation factor Xa, or urokinase-type plasminogen activator., similarity: Belongs to the serpin family., similarity: Belongs to the serpin family. Ov-serpin subfamily., tissue specificity: Expressed in many tissues, including brain, bone marrow, lymph node, heart, lung, liver, pancreas, testis, ovary, and


intestine.,




+86-27-59760950 ELKbio@ELKbiotech.com

www.elkbiotech.com






Western Blot analysis of various cells using Serpin B12 Polyclonal Antibody




Western blot analysis of SERPINB12 Antibody. The lane on the right is blocked with the SERPINB12 peptide.



Western blot analysis of the lysates from HeLa cells using SERPINB12 antibody.







Immunohistochemical analysis of paraffin-embedded human cervical carcinoma. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).



+86-27-59760950