

## Skp1 p19 rabbit pAb

### Cat No.:ES3449

For research use only

#### Overview

Product Name	Skp1 p19 rabbit pAb
Host species	Rabbit
Applications	WB;IHC;IF;ELISA
Species Cross-Reactivity	Human;Mouse;Rat;Monkey
<b>Recommended dilutions</b>	Western Blot: 1/500 - 1/2000. Immunohistochemistry:
	1/100 - 1/300. ELISA: 1/5000. Not yet tested in other
	applications.
Immunogen	The antiserum was produced against synthesized
	peptide derived from human SKP1A/p19. AA
	range:41-90
Specificity	Skp1 p19 Polyclonal Antibody detects endogenous
	levels of Skp1 p19 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	S-phase kinase-associated protein 1
Gene Name	SKP1
Cellular localization	nucleus, nucle oplasm, cytoplasm, centrosome, cytosol, SCF
	ubiquitin ligase complex,Cul7-RING ubiquitin ligase
	complex,PcG protein complex,extracellular exosome,
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	19kD
Human Gene ID	6500
Human Swiss-Prot	P63208
Number	
Alternative Names	SKP1; EMC19; OCP2; SKP1A; TCEB1L; S-phase
	kinase-associated protein 1; Cyclin-A/CDK2-associated
	protein p19; Organ of Corti protein 2; OCP-2; Organ of
	Corti protein II; OCP-II; RNA polymerase II elongation



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com


23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C




Background


#### factor-like protein; SIII; Transcr

This gene encodes a component of SCF complexes, which are composed of this protein, cullin 1, a ring-box protein, and one member of the F-box family of proteins. This protein binds directly to the F-box motif found in F-box proteins. SCF complexes are involved in the regulated ubiquitination of specific protein substrates, which targets them for degradation by the proteosome. Specific F-box proteins recognize different target protein(s), and many specific SCF substrates have been identified including regulators of cell cycle progression and development. Studies have also characterized the protein as an RNA polymerase II elongation factor. Alternative splicing of this gene results in two transcript variants. A related pseudogene has been identified on chromosome 7. [provided by RefSeq, Jul 2008],





Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using SKP1A/p19 Antibody. The picture on the right is blocked with the synthesized peptide.





+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



Western blot analysis of lysates from COS7 cells, using SKP1A/p19 Antibody. The lane on the right is blocked with the synthesized peptide.

# SKP1A/p19-



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C