

Trypsin-3 rabbit pAb

Cat No.:ES3649

For research use only

Overview

Product Name	Trypsin-3 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not
	yet tested in other applications.
Immunogen	The antiserum was produced against synthesized
	peptide derived from human Trypsin-3. AA
	range:117-166
Specificity	Trypsin-3 Polyclonal Antibody detects endogenous
	levels of Trypsin-3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Trypsin-3
Gene Name	PRSS3
Cellular localization	Secreted .
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	32kD
Human Gene ID	5646
Human Swiss-Prot Number	P35030
Alternative Names	PRSS3; PRSS4; TRY3; TRY4; Trypsin-3; Brain
	trypsinogen; Mesotrypsinogen; Serine protease 3;
	Serine protease 4; Trypsin III; Trypsin IV
Background	This gene encodes a trypsinogen, which is a member
	of the trypsin family of serine proteases. This
	enzyme is expressed in the brain and pancreas and
	is resistant to common trypsin inhibitors. It is active
	on peptide linkages involving the carboxyl group of



+86-27-59760950

ELKbio@ELKbiotech.com


www.elkbiotech.com


23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C




lysine or arginine. This gene is localized to the locus of T cell receptor beta variable orphans on chromosome 9. Four transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Oct 2010],

Western Blot analysis of various cells using Trypsin-3 Polyclonal Antibody. Secondary antibody(catalog#:RS0002) was diluted at 1:20000









+86-27-59760950

m www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C