

CD1D rabbit pAb

Cat No.: ES3969

For research use only

Overview

Product Name CD1D rabbit pAb

Host species Rabbit

Applications WB;IHC;IF;ELISA Species Cross-Reactivity Human;Rat;Mouse;

Recommended dilutions Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300

ELISA: 1/20000. Not yet tested in other applications.

Immunogen Synthesized peptide derived from

Antigen-presenting glycoprotein CD1d at AA range:

161-210

Specificity CD1D Polyclonal Antibody detects endogenous

levels of CD1D protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Antigen-presenting glycoprotein CD1d

Gene Name CD1D

Cell ular localization Cell membrane ; Single-pass type I membrane

protein . Basolateral cell membrane ; Single-pass type I membrane protein . Endosome membrane ; Single-pass type I membrane protein . Lysosome membrane ; Single-pass type I membrane protein . Endoplasmic reticulum membrane ; Single-pass type

I membrane protein . Subject to intracellular

trafficking between the cell membrane, endosomes

and lysosomes...

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 37kD
Human Gene ID 912
Human Swiss-Prot Number P15813



+86-27-59760950 ELKbio@ELKbiotech.com www.elkbiotech.com




Alternative Names

Background

CD1D; Antigen-presenting glycoprotein CD1d; R3G1; CD1d

This gene encodes a divergent member of the CD1 family of transmembrane glycoproteins, which are structurally related to the major histocompatibility complex (MHC) proteins and form heterodimers with beta-2-microglobulin. The CD1 proteins mediate the presentation of primarily lipid and glycolipid antigens of self or microbial origin to T cells. The human genome contains five CD1 family genes organized in a cluster on chromosome 1. The CD1 family members are thought to differ in their cellular localization and specificity for particular lipid ligands. The protein encoded by this gene localizes to late endosomes and lysosomes via a tyrosine-based motif in the cytoplasmic tail. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2016],



Western Blot analysis of K562 cells using CD1D Polyclonal Antibody. Antibody was diluted at 1:500. Secondary antibody(catalog#:RS0002) was diluted at 1:20000

Western blot analysis of SH-SY5Y lysis using CD1D antibody. Antibody was diluted at 1:500. Secondary antibody(catalog#:RS0002) was diluted at 1:20000



+86-27-59760950

CD1D

35---25---


ELKbio@ELKbiotech.com

www.elkbiotech.com





Immunohistochemical analysis of paraffin-embedded human-lymph, antibody was diluted at 1:100



+86-27-59760950

Immunohistochemical analysis of paraffin-embedded human-lymph, antibody was diluted at 1:100

