

## CAF-1 p150 rabbit pAb

**Cat No.:ES4412** 

For research use only

## Overview

Product Name CAF-1 p150 rabbit pAb

Host species Rabbit
Applications WB;ELISA
Species Cross-Reactivity Human;Mouse

Recommended dilutions Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not

yet tested in other applications.

Immunogen Synthesized peptide derived from CAF-1 p150 . at

AA range: 300-380

**Specificity** CAF-1 p150 Polyclonal Antibody detects endogenous

levels of CAF-1 p150 protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

**Storage** Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Chromatin assembly factor 1 subunit A

Gene Name CHAF1A

Cellular localization Nucleus . DNA replication foci.

**Purification** The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 150kD
Human Gene ID 10036
Human Swiss-Prot Number Q13111

Alternative Names CHAF1A; CAF; CAF1P150; Chromatin assembly factor

1 subunit A; CAF-1 subunit A; Chromatin assembly factor I p150 subunit; CAF-I 150 kDa subunit; CAF-I

p150; hp150

**Background** Chromatin assembly factor I (CAF1) is a nuclear


complex consisting of p50, p60 (CHAF1B; MIM 601245), and p150 (CHAF1A) subunits that

assembles histone octamers onto replicating DNA in

vitro (Kaufman et al., 1995 [PubMed







## 7600578]).[supplied by OMIM, Mar 2008],

Western Blot analysis of K562 cells using CAF-1 p150 Polyclonal Antibody. Secondary antibody(catalog#:RS0002) was diluted at 1:20000 cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003,Inventbiotech,MN,USA).



+86-27-59760950