

Hexim1 rabbit pAb

Cat No.: ES4522

For research use only

Overview

Product Name Hexim1 rabbit pAb

Host species Rabbit

Applications WB;IHC;IF;ELISA Species Cross-Reactivity Human;Mouse;Rat

Recommended dilutions Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA:

1/10000. Not yet tested in other applications.

The antiserum was produced against synthesized

Immunogen The antiserum was produced against synthesize peptide derived from human HEXIM1. AA

range:181-230

Specificity Hexim1 Polyclonal Antibody detects endogenous

levels of Hexim1 protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Protein HEXIM1

Gene Name HEXIM1

Cellular localization Nucleus . Cytoplasm . Binds alpha-importin and is

mostly nuclear (PubMed:16362050).

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 40kD
Human Gene ID 10614
Human Swiss-Prot Number 094992

Alternative Names HEXIM1; CLP1; EDG1; HIS1; MAQ1; Protein HEXIM1;

Cardiac lineage protein 1; Estrogen down-regulated

gene 1 protein; Hexamethylene

bis-acetamide-inducible protein 1; Menage a quatre


protein 1



+86-27-59760950 ELKbio@ELKbiotech.com www.elkbiotech.com




Background





Expression of this gene is induced by hexamethylene-bis-acetamide in vascular smooth muscle cells. This gene has no introns. [provided by RefSeq, Jul 2008],

Western Blot analysis of RAT-HEART cells using Hexim1 Polyclonal Antibody diluted at 1:2000

Immunofluorescence analysis of HeLa cells, using HEXIM1 Antibody. The picture on the right is blocked with the synthesized peptide.








Western blot analysis of lysates from NIH/3T3 cells, using HEXIM1 Antibody. The lane on the right is blocked with the synthesized peptide.







Western blot analysis of the lysates from HeLa cells using $\ensuremath{\mathsf{HEXIM1}}$ antibody.

