

Ataxin-2L rabbit pAb

Cat No.: ES4628

For research use only

Overview

Product Name Ataxin-2L rabbit pAb

Host species Rabbit
Applications WB;ELISA

Species Cross-Reactivity Human;Rat;Mouse;

Recommended dilutions Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not

yet tested in other applications.

Immunogen The antiserum was produced against synthesized

peptide derived from human ATXN2L. AA

range:571-620

Specificity Ataxin-2L Polyclonal Antibody detects endogenous

levels of Ataxin-2L protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Ataxin-2-like protein

Gene Name ATXN2L

Cellular localization Membrane ; Peripheral membrane protein .

Cytoplasm . Nucleus speckle . Cytoplasmic granule . Predominantly cytoplasmic but is also detected in nuclear speckles (PubMed:23209657). Component of cytoplasmic stress granules (PubMed:23209657).

Inhibition of me

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 113kD
Human Gene ID 11273

Human Swiss-Prot Number Q8WWM7

Alternative Names ATXN2L; A2D; A2LG; A2LP; A2RP; Ataxin-2-like

protein; Ataxin-2 domain protein; Ataxin-2-related

protein




+86-27-59760950 ELKbio@ELKbiotech.com

www.elkbiotech.com




Background


This gene encodes an ataxin type 2 related protein of unknown function. This protein is a member of the spinocerebellar ataxia (SCAs) family, which is associated with a complex group of neurodegenerative disorders. Several alternatively spliced transcripts encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008],



Western Blot analysis of 293T cells using Ataxin-2L Polyclonal Antibody diluted at 1:2000



Western blot analysis of lysates from NIH/3T3 and RAW264.7 cells, using ATXN2L Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from 293 cells using ATXN2L antibody.



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com