



ABHD12B rabbit pAb

Cat No.:ES4860

For research use only

Overview


Product Name	ABHD12B rabbit pAb
Host species	Rabbit
Applications	IHC;IF;ELISA
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Immunogen	The antiserum was produced against synthesized peptide derived from human ABHD12B. AA range:281-330
Specificity	ABHD12B Polyclonal Antibody detects endogenous levels of ABHD12B protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Abhydrolase domain-containing protein 12B
Gene Name	ABHD12B
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	40kD
Human Gene ID	145447
Human Swiss-Prot Number	Q7Z5M8
Alternative Names	ABHD12B; C14orf29; Abhydrolase domain-containing protein 12B
Background	similarity:Belongs to the serine esterase family.,





ELK Biotechnology

Immunohistochemistry analysis of paraffin-embedded human colon carcinoma, using ABHD12B Antibody. The picture on the right is blocked with the synthesized peptide.



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei, P.R.C