

PDIK1L rabbit pAb

Cat No.: ES4886

For research use only

Overview

Product Name PDIK1L rabbit pAb

Host species Rabbit

ApplicationsWB;IHC;IF;ELISASpecies Cross-ReactivityHuman;Mouse

Recommended dilutions Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.

Immunogen The antiserum was produced against synthesized

peptide derived from human PDIK1L. AA

range:101-150

Specificity PDIK1L Polyclonal Antibody detects endogenous

levels of PDIK1L protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Serine/threonine-protein kinase PDIK1L

Gene Name PDIK1L **Cellular localization** Nucleus .

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

ClonalityPolyclonalConcentration1 mg/mlObserved band39kDHuman Gene ID149420Human Swiss-Prot NumberQ8N165

Alternative Names PDIK1L; CLIK1L; Serine/threonine-protein kinase

PDIK1L; PDLIM1-interacting kinase 1-like

Background catalytic activity:ATP + a protein = ADP + a

phosphoprotein., similarity: Belongs to the protein

kinase superfamily. Ser/Thr protein kinase family., similarity: Contains 1 protein kinase

domain.,tissue specificity:Expressed in liver, kidney,



+86-27-59760950 ELKbio@ELKbiotech.com

www.elkbiotech.com




K562 138--100--70--55--40-- ____ PDIK1L

35---25---


15---

pancreas, spleen, thymus and prostate.,

Western Blot analysis of K562 cells using PDIK1L Polyclonal Antibody cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003,Inventbiotech,MN,USA).



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using PDIK1L Antibody. The picture on the right is blocked with the synthesized peptide.



+86-27-59760950

Western blot analysis of lysates from K562 cells, using PDIK1L Antibody. The lane on the right is blocked with the synthesized peptide.



23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei, P.R.C