

CYP2B6 rabbit pAb

Cat No.:ES4937

For research use only

Overview

Product Name	CYP2B6 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human; Rat; Mouse;
Recommended dilutions	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not
	yet tested in other applications.
Immunogen	Synthesized peptide derived from the Internal
	region of human CYP2B6.
Specificity	CYP2B6 Polyclonal Antibody detects endogenous
	levels of CYP2B6 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Cytochrome P450 2B6
Gene Name	CYP2B6
Cellular localization	Endoplasmic reticulum membrane; Peripheral
	membrane protein. Microsome membrane;
	Peripheral membrane protein.
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	55kD
Human Gene ID	1565
Human Swiss-Prot Number	P20813
Alternative Names	CYP2B6; Cytochrome P450 2B6; 1; 4-cineole
	2-exo-monooxygenase; CYPIIB6; Cytochrome P450
	IIB1
Background	This gene, CYP2B6, encodes a member of the
	cytochrome P450 superfamily of enzymes. The
	cytochrome P450 proteins are monooxygenases
	which catalyze many reactions involved in drug
_	



+86-27-59760950


ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



metabolism and synthesis of cholesterol, steroids and other lipids. This protein localizes to the endoplasmic reticulum and its expression is induced by phenobarbital. The enzyme is known to metabolize some xenobiotics, such as the anti-cancer drugs cyclophosphamide and ifosphamide. Transcript variants for this gene have been described; however, it has not been resolved whether these transcripts are in fact produced by this gene or by a closely related pseudogene, CYP2B7. Both the gene and the pseudogene are located in the middle of a CYP2A pseudogene found in a large cluster of cytochrome P450 genes from the CYP2A, CYP2B and CYP2F subfamilies on chromosome 19q. [provided by RefSeq, Jul 2008],



Western Blot analysis of 293 cells using CYP2B6 Polyclonal Antibody



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C