

ZNF600 rabbit pAb

Cat No.:ES4971

For research use only

Overview

Product Name	ZNF600 rabbit pAb	
Host species	Rabbit	
Applications	IHC;IF;ELISA	
Species Cross-Reactivity	Human;Rat;Mouse;	
Recommended dilutions	Immunohistochemistry: 1/100 - 1/300. ELISA:	
	1/10000. Not yet tested in other applications.	
Immunogen	The antiserum was produced against synthesized	
	peptide derived from human ZNF600. AA	
	range:250-300	
Specificity	ZNF600 Polyclonal Antibody detects endogenous	
	levels of ZNF600 protein.	
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and	
	0.02% sodium azide.	
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.	
Protein Name	Zinc finger protein 600	
Gene Name	ZNF600	
Cellular localization	Nucleus .	
Purification	The antibody was affinity-purified from rabbit	
	antiserum by affinity-chromatography using	
	epitope-specific immunogen.	
Clonality	Polyclonal	
Concentration	1 mg/ml	
Observed band		
Human Gene ID	162966	
Human Swiss-Prot Number	Q6ZNG1	
Alternative Names	ZNF600; Zinc finger protein 600	
Background	function:May be involved in transcriptional	
	regulation.,PTM:Phosphorylated upon DNA damage,	
	probably by ATM or ATR., similarity:Belongs to the	
	krueppel C2H2-type zinc-finger protein	
	family., similarity: Contains 20 C2H2-type zinc	
	fingers.,	




+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C





Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C