

## ZNF596 rabbit pAb

Cat No.: ES4994

For research use only

## Overview

Product Name ZNF596 rabbit pAb

Host species Rabbit

Applications WB;IHC;IF;ELISA Species Cross-Reactivity Human;Rat;Mouse;

**Recommended dilutions** Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.

Immunogen The antiserum was produced against synthesized

peptide derived from human ZNF596. AA

range:261-310

**Specificity** ZNF596 Polyclonal Antibody detects endogenous

levels of ZNF596 protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

**Storage** Store at -20°C. Avoid repeated freeze-thaw cycles.

**Protein Name** Zinc finger protein 596

**Gene Name** ZNF596 **Cellular localization** Nucleus .

**Purification** The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

ClonalityPolyclonalConcentration1 mg/mlObserved band58kDHuman Gene ID169270Human Swiss-Prot NumberQ8TC21

Alternative Names ZNF596; Zinc finger protein 596

**Background** function: May be involved in transcriptional

regulation., similarity: Belongs to the krueppel

C2H2-type zinc-finger protein family., similarity: Contains 1 KRAB

domain., similarity: Contains 11 C2H2-type zinc


fingers.,



+86-27-59760950 ELKbio@ELKbiotech.com


www.elkbiotech.com





Western Blot analysis of Jurkat cells using ZNF596
Polyclonal Antibody. Secondary
antibody(catalog#:RS0002) was diluted at 1:20000 cells
nucleus extracted by Minute TM Cytoplasmic and Nuclear
Fractionation kit (SC-003,Inventbiotech,MN,USA).

Immunohistochemistry analysis of paraffin-embedded human brain tissue, using ZNF596 Antibody. The picture on the right is blocked with the synthesized peptide.



ZNF596 -- -- 48 -- 34 -- 26 -- 19 (kD)

Western blot analysis of lysates from Jurkat cells, using ZNF596 Antibody. The lane on the right is blocked with the synthesized peptide.

