

## Eps15 (phospho Tyr849) rabbit pAb

Cat No.: ES5140

For research use only

## Overview

Product Name Eps15 (phospho Tyr849) rabbit pAb

Host species Rabbit
Applications WB;ELISA

Species Cross-Reactivity Human; Mouse; Rat

Recommended dilutions Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not

yet tested in other applications.

Immunogen Synthesized phospho-peptide around the

phosphorylation site of human Eps15 (phospho

Tyr849)

**Specificity** Phospho-Eps15 (Y849) Polyclonal Antibody detects

endogenous levels of Eps15 protein only when

phosphorylated at Y849.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Epidermal growth factor receptor substrate 15

Gene Name EPS15

**Cellular localization** Cytoplasm. Cell membrane; Peripheral membrane

protein; Cytoplasmic side. Membrane, clathrin-coated pit. Recruited to the plasma membrane upon EGFR activation and localizes to

coated pits. Colocalizes with UBQLN1 in ubiquitin-rich cytoplasmic aggregates th

**Purification** The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 100kD
Human Gene ID 2060
Human Swiss-Prot Number P42566

Alternative Names EPS15; AF1P; Epidermal growth factor receptor

substrate 15; Protein Eps15; Protein AF-1p




+86-27-59760950 ELKbio@ELKbiotech.com

www.elkbiotech.com



**Background** 



This gene encodes a protein that is part of the EGFR pathway. The protein is present at clatherin-coated pits and is involved in receptor-mediated endocytosis of EGF. Notably, this gene is rearranged with the HRX/ALL/MLL gene in acute myelogeneous leukemias. Alternative splicing results in multiple transcript variants encoding distinct isoforms. [provided by RefSeq, May 2009],

Western Blot analysis of JK cells using Phospho-Eps15 (Y849) Polyclonal Antibody



+86-27-59760950