

## ERF (phospho Thr526) rabbit pAb

Cat No.: ES5168

For research use only

## Overview

**Immunogen** 

Product Name ERF (phospho Thr526) rabbit pAb

Host species Rabbit

**Applications** WB;IHC;IF;ELISA **Species Cross-Reactivity** Human;Mouse

**Recommended dilutions** Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications. The antiserum was produced against synthesized

peptide derived from human ERF around the phosphorylation site of Thr526. AA range:492-541

**Specificity** Phospho-ERF (T526) Polyclonal Antibody detects

endogenous levels of ERF protein only when

phosphorylated at T526.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name ETS domain-containing transcription factor ERF

Gene Name ERF
Cellular localization Nucleus.

**Purification** The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 58kD
Human Gene ID 2077
Human Swiss-Prot Number P50548

Alternative Names ERF; ETS domain-containing transcription factor ERF;

Ets2 repressor factor; PE-2


**Background** ETS2 is a transcription factor and protooncogene

involved in development, apoptosis, and the regulation of telomerase. The protein encoded by







this gene binds to the ETS2 promoter and is a strong repressor of ETS2 transcription. Several transcript variants encoding two different isoforms have been found for this gene. [provided by RefSeq, Aug 2015],



Enzyme-Linked Immunosorbent Assay (Phospho-ELISA) for Immunogen Phosphopeptide (Phospho-left) and Non-Phosphopeptide (Phospho-right), using ERF (Phospho-Thr526) Antibody



Immunofluorescence analysis of HeLa cells, using ERF (Phospho-Thr526) Antibody. The picture on the right is blocked with the phospho peptide.



+86-27-59760950

Immunohistochemistry analysis of paraffin-embedded human lung carcinoma, using ERF (Phospho-Thr526) Antibody. The picture on the right is blocked with the phospho peptide.

