

FoxJ3 rabbit pAb

Cat No.: ES5268

For research use only

Overview

Product Name FoxJ3 rabbit pAb

Host species Rabbit

Applications IHC;IF;ELISA **Species Cross-Reactivity** Human;Mouse

Recommended dilutions Immunohistochemistry: 1/100 - 1/300. ELISA:

1/20000. Not yet tested in other applications.

Immunogen The antiserum was produced against synthesized

peptide derived from human FOXJ3. AA

range:191-240

Specificity FoxJ3 Polyclonal Antibody detects endogenous

levels of FoxJ3 protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Forkhead box protein J3

Gene Name FOXJ3 **Cellular localization** Nucleus .

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal Concentration 1 mg/ml

Observed band

Human Gene ID 22887 Human Swiss-Prot Number Q9UPW0


+86-27-59760950

Alternative Names FOXJ3; KIAA1041; Forkhead box protein J3
Background similarity:Contains 1 fork-head DNA-binding

domain.,







Immunohistochemistry analysis of paraffin-embedded human brain tissue, using FOXJ3 Antibody. The picture on the right is blocked with the synthesized peptide.



+86-27-59760950