

Rabphilin-3A (phospho Ser237) rabbit pAb

Cat No.:ES5269

For research use only

Overview

Immunogen

Purification

Product Name Rabphilin-3A (phospho Ser237) rabbit pAb

Host species Rabbit

Applications WB;IHC;IF;ELISA **Species Cross-Reactivity** Human;Mouse;Rat

Recommended dilutions Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications. The antiserum was produced against synthesized

peptide derived from human Rabphilin 3A around

the phosphorylation site of Ser237. AA

range:203-252

Specificity Phospho-Rabphilin-3A (S237) Polyclonal Antibody

detects endogenous levels of Rabphilin-3A protein

only when phosphorylated at S237.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Rabphilin-3A
Gene Name RPH3A

Cellular localization Cytoplasmic vesicle, secretory vesicle, synaptic

vesicle membrane . Cell projection, dendritic spine . Cell junction, synapse, postsynaptic cell membrane .

Membrane; Peripheral membrane protein.

The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 75kD
Human Gene ID 22895
Human Swiss-Prot Number Q9Y2J0


Alternative Names RPH3A; KIAA0985; Rabphilin-3A; Exophilin-1






Background


Exocytosis of neurotransmitters and hormones is fundamental to synaptic neurotransmission and cell-cell communication. RAB3A (MIM 179390) is a small G protein that is thought to act at late stages of exocytosis, and RPH3A is a RAB3A effector (Lin et al., 2007 [PubMed 17149709]).[supplied by OMIM, Jul 2008],



Enzyme-Linked Immunosorbent Assay (Phospho-ELISA) for Immunogen Phosphopeptide (Phospho-left) and Non-Phosphopeptide (Phospho-right), using Rabphilin 3A (Phospho-Ser237) Antibody



Immunofluorescence analysis of HeLa cells, using Rabphilin 3A (Phospho-Ser237) Antibody. The picture on the right is blocked with the phospho peptide.




+86-27-59760950

Immunohistochemistry analysis of paraffin-embedded human brain, using Rabphilin 3A (Phospho-Ser237) Antibody. The picture on the right is blocked with the phospho peptide.







Western blot analysis of lysates from RAW264.7 cells treated with Calyculin 100nM 15', using Rabphilin 3A (Phospho-Ser237) Antibody. The lane on the right is blocked with the phospho peptide.



+86-27-59760950