

## FoxE3 rabbit pAb

## Cat No.:ES5288

For research use only

## Overview

Product Name	FoxE3 rabbit pAb
Host species	Rabbit
Applications	WB;IHC;IF;ELISA
Species Cross-Reactivity	Human;Mouse;Rat
<b>Recommended dilutions</b>	Western Blot: 1/500 - 1/2000.
	Immunohistochemistry: 1/100 - 1/300. ELISA:
	1/20000. Not yet tested in other applications.
Immunogen	The antiserum was produced against synthesized
	peptide derived from human FOXE3. AA
	range:81-130
Specificity	FoxE3 Polyclonal Antibody detects endogenous
	levels of FoxE3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Forkhead box protein E3
Gene Name	FOXE3
<b>Cellular localization</b>	Nucleus .
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	33kD
Human Gene ID	2301
Human Swiss-Prot Number	Q13461
Alternative Names	FOXE3; FKHL12; FREAC8; Forkhead box protein E3;
	Forkhead-related protein FKHL12; Forkhead-related
	transcription factor 8; FREAC-8
Background	This intronless gene belongs to the forkhead family
	of transcription factors, which is characterized by a
	distinct forkhead domain. The protein encoded
	functions as a lens-specific transcription factor and



+86-27-59760950

ELKbio@ELKbiotech.com


www.elkbiotech.com


23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



plays an important role in vertebrate lens formation. Mutations in this gene are associated with anterior segment mesenchymal dysgenesis and congenital primary aphakia. [provided by RefSeq, Dec 2009],

Western Blot analysis of K562 cells using FoxE3 Polyclonal Antibody cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003,Inventbiotech,MN,USA).





K562

FoxE3

138-

70--

55---40---35---

25-

15----

K562 -- 117 -- 85 FOXE3---- 48 -- 34 -- 26 -- 19 (kD)

Western blot analysis of lysates from K562 cells, using FOXE3 Antibody. The lane on the right is blocked with the synthesized peptide.



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C