

Palladin rabbit pAb

Cat No.: ES5291

For research use only

Overview

Product Name Palladin rabbit pAb

Host species Rabbit
Applications IHC;IF;ELISA
Species Cross-Reactivity Human;Mouse

Recommended dilutions Immunohistochemistry: 1/100 - 1/300. ELISA:

1/20000. Not yet tested in other applications.

Immunogen Synthesized peptide derived from Palladin . at AA

range: 450-530

Specificity Palladin Polyclonal Antibody detects endogenous

levels of Palladin protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Palladin
Gene Name PALLD

Cellular localization Cytoplasm, cytoskeleton . Cell junction, focal

adhesion . Cytoplasm, myofibril, sarcomere, Z line . Cell projection, ruffle . Cell projection, podosome . Cell projection, lamellipodium . Cell projection, axon . Cell projection, growth cone . Localizes to The antibody was affinity-purified from rabbit

Purification The antibody was affinity-purified from rabbi antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 150kD
Human Gene ID 23022
Human Swiss-Prot Number Q8WX93

Alternative Names PALLD; KIAA0992; CGI-151; Palladin; SIH002;

Sarcoma antigen NY-SAR-77

Background This gene encodes a cytoskeletal protein that is

required for organizing the actin cytoskeleton. The

protein is a component of actin-containing



+86-27-59760950 ELKbio@ELKbiotech.com

www.elkbiotech.com




microfilaments, and it is involved in the control of cell shape, adhesion, and contraction.

Polymorphisms in this gene are associated with a susceptibility to pancreatic cancer type 1, and also with a risk for myocardial infarction. Alternative splicing results in multiple transcript variants.

[provided by RefSeq, Oct 2009],

Immunohistochemical analysis of paraffin-embedded human Colon cancer. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).



+86-27-59760950

