

## Serpin A11 rabbit pAb

Cat No.: ES5433

For research use only

## Overview

Product Name Serpin A11 rabbit pAb

Host species Rabbit
Applications WB;IHC

**Species Cross-Reactivity** Human; Rat; Mouse;

Recommended dilutions WB 1:500-2000;IHC-p 1:50-300

Immunogen The antiserum was produced against synthesized

peptide derived from human SERPINA11. AA

range:191-240

**Specificity** Serpin A11 Polyclonal Antibody detects endogenous

levels of Serpin A11 protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

**Storage** Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Serpin A11
Gene Name SERPINA11
Cellular localization Secreted.

**Purification** The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 47kD
Human Gene ID 256394
Human Swiss-Prot Number Q86U17


+86-27-59760950

Alternative Names SERPINA11; Serpin A11


**Background** similarity:Belongs to the serpin family.,







Western blot analysis of lysates from HT-29 cells, using SERPINA11 Antibody. The lane on the right is blocked with the synthesized peptide.



Immunohistochemical analysis of paraffin-embedded human liver cancer. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).

