

## SRPK3 rabbit pAb

## Cat No.:ES5502

For research use only

## Overview

Product Name	SRPK3 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Mouse;Rat
<b>Recommended dilutions</b>	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not
	yet tested in other applications.
Immunogen	The antiserum was produced against synthesized
-	peptide derived from human SRPK3. AA
	range:261-310
Specificity	SRPK3 Polyclonal Antibody detects endogenous
. ,	levels of SRPK3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	SRSF protein kinase 3
Gene Name	SRPK3
Cellular localization	
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	60kD
Human Gene ID	26576
Human Swiss-Prot Number	Q9UPE1
Alternative Names	SRPK3; MSSK1; STK23; SRSF protein kinase 3;
	Muscle-specific serine kinase 1; MSSK-1;
	Serine/arginine-rich protein-specific kinase 3;
	SR-protein-specific kinase 3;
	Serine/threonine-protein kinase 23
Background	This gene encodes a protein kinase similar to a
	protein kinase which is specific for the SR
	(serine/arginine-rich domain) family of splicing
_	, , , , , , , , , , , , , , , , , , ,



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



293 HeLa

SRPK3

(kD)

117-

85-

48-

34-

26-

19-

--150 --100 --75

--50 --37

--25 --20

--15 (kd)

SRPK3

factors. A highly similar protein has been shown to play a role in muscle development in mice. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2009],

Western blot analysis of SRPK3 Antibody. The lane on the right is blocked with the SRPK3 peptide.

Western blot analysis of the lysates from HeLa cells using SRPK3 antibody.



ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C