

Cleaved-GGT4 HC (T472) rabbit pAb

Cat No.:ES5516

For research use only

Overview

Product Name Cleaved-GGT4 HC (T472) rabbit pAb

Host species Rabbit
Applications WB;ELISA

Species Cross-Reactivity Human; Mouse; Rat

Recommended dilutions Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not

yet tested in other applications.

Immunogen The antiserum was produced against synthesized

peptide derived from human

Gamma-glutamyltransferase 4. AA range:423-472

Specificity Cleaved-GGT4 HC (T472) Polyclonal Antibody

detects endogenous levels of fragment of activated GGT4 HC protein resulting from cleavage adjacent to

T472.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Gamma-glutamyltransferase 7

Gene Name GGT7

Cellular localization Membrane; Single-pass type II membrane protein.

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 50kD
Human Gene ID 2686
Human Swiss-Prot Number O9UJ14

Alternative Names GGT7; GGTL3; GGTL5; Gamma-glutamyltransferase

7; GGT 7; Gamma-glutamyltransferase-like 3;

Gamma-glutamyltransferase-like 5;

Gamma-glutamyltranspeptidase 7; Glutathione

hydrolase 7

Background This gene is a member of a gene family that encodes



+86-27-59760950 ELKbio@ELKbiotech.com www.elkbiotech.com



enzymes involved in both the metabolism of glutathione and in the transpeptidation of amino acids. Changes in the activity of gamma-glutamyltransferase may signal preneoplastic or toxic conditions in the liver or kidney. The protein encoded by this gene consists of a heavy and a light chain, and it can interact with CT120, a plasma membrane-associated protein that is possibly involved in lung carcinogenesis. [provided by RefSeq, Jul 2008],

Western blot analysis of lysates from Jurkat cells, treated with etoposide 25uM 24h, using Gamma-glutamyltransferase 4 (heavy chain,Cleaved-Thr472) Antibody. The lane on the right is blocked with the synthesized peptide.

γ-glutamyltransferase 4 (heavy chain,Cleaved-Thr472

