

## CYP4V2 rabbit pAb

Cat No.: ES5608

For research use only

## Overview

Product Name CYP4V2 rabbit pAb

Host species Rabbit

Applications WB;IHC;IF;ELISA Species Cross-Reactivity Human;Rat;Mouse;

**Recommended dilutions** Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.

**Immunogen** Synthesized peptide derived from the Internal

region of human CYP4V2.

**Specificity** CYP4V2 Polyclonal Antibody detects endogenous

levels of CYP4V2 protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

**Storage** Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Cytochrome P450 4V2

Gene Name CYP4V2

**Cellular localization** Endoplasmic reticulum membrane ; Single-pass

membrane protein.

**Purification** The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

ClonalityPolyclonalConcentration1 mg/mlObserved band60kDHuman Gene ID285440Human Swiss-Prot NumberQ6ZWL3

Alternative Names CYP4V2; Cytochrome P450 4V2

**Background** This gene encodes a member of the cytochrome

P450 hemethiolate protein superfamily which are involved in oxidizing various substrates in the metabolic pathway. It is implicated in the metabolism of fatty acid precursors into n-3

polyunsaturated fatty acids. Mutations in this gene




+86-27-59760950 ELKbio@ELKbiotech.com

www.elkbiotech.com



result in Bietti crystalline corneoretinal dystrophy. [provided by RefSeq, Jul 2008],



Immunohistochemical analysis of paraffin-embedded human Gastric adenocarcinoma. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).



+86-27-59760950