

GluR4 (phospho Ser862) rabbit pAb

Cat No.: ES5637

For research use only

Overview

Product Name GluR4 (phospho Ser862) rabbit pAb

Host species Rabbit

Applications WB;IHC;IF;ELISA Species Cross-Reactivity Human;Mouse;Rat

Recommended dilutions Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.

Immunogen The antiserum was produced against synthesized

peptide derived from human GluR4 around the phosphorylation site of Ser862. AA range:828-877

Specificity Phospho-GluR4 (S862) Polyclonal Antibody detects

endogenous levels of GluR4 protein only when

phosphorylated at S862.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Glutamate receptor 4

Gene Name GRIA4

Cellular localization Cell membrane; Multi-pass membrane protein. Cell

junction, synapse, postsynaptic cell membrane; Multi-pass membrane protein. Cell projection, dendrite. Interaction with CNIH2, CNIH3 and PRKCG

promotes cell surface expression. .

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 100kD
Human Gene ID 2893
Human Swiss-Prot Number P48058

Alternative Names GRIA4; GLUR4; Glutamate receptor 4; GluR-4;

GluR4; AMPA-selective glutamate receptor 4;



+86-27-59760950 ELKbio@ELKbiotech.com

www.elkbiotech.com




Background

GluR-D; Glutamate receptor ionotropic; AMPA 4; GluA4


Glutamate receptors are the predominant excitatory neurotransmitter receptors in the mammalian brain and are activated in a variety of normal neurophysiologic processes. These receptors are heteromeric protein complexes composed of multiple subunits, arranged to form ligand-gated ion channels. The classification of glutamate receptors is based on their activation by different pharmacologic agonists. The subunit encoded by this gene belongs to a family of AMPA

(alpha-amino-3-hydroxy-5-methyl-4-isoxazole propionate)-sensitive glutamate receptors, and is subject to RNA editing (AGA->GGA; R->G). Alternative splicing of this gene results in transcript variants encoding different isoforms, which may vary in their signal transduction properties. Some haplotypes of this gene show a positive association with schizophrenia. [provided by RefSeq, Jul 2008],

Immunohistochemistry analysis of paraffin-embedded human brain, using GluR4 (Phospho-Ser862) Antibody. The picture on the right is blocked with the phospho peptide.



Western blot analysis of lysates from HepG2 cells treated with Forskolin 40nM 30', using GluR4 (Phospho-Ser862) Antibody. The lane on the right is blocked with the phospho peptide.



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com







+86-27-59760950