

## GRF-1 (phospho Tyr1087) rabbit pAb

## Cat No.:ES5668

For research use only

## Overview

Product Name	GRF-1 (phospho Tyr1087) rabbit pAb	
Host species	Rabbit	
Applications	WB;ELISA	
Species Cross-Reactivity	Human;Mouse;Rat	
Recommended dilutions	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not	
	yet tested in other applications.	
Immunogen	Synthesized phospho-peptide around the	
	phosphorylation site of human GRF-1 (phospho	
	Tyr1087)	
Specificity	Phospho-GRF-1 (Y1087) Polyclonal Antibody detects	
	endogenous levels of GRF-1 protein only when	
	phosphorylated at Y1087.	
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and	
	0.02% sodium azide.	
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.	
Protein Name	Rho GTPase-activating protein 35	
Gene Name	ARHGAP35	
Cellular localization	Cytoplasm, cytoskeleton, cilium basal body .	
	Cytoplasm . Nucleus . Cell membrane . In response	
	to integrins and SDC4 and upon phosphorylation by	
	PKC, relocalizes from the cytoplasm to regions of	
	plasma membrane ruffling where it colocalizes with	
	polymeriz	
Purification	The antibody was affinity-purified from rabbit	
	antiserum by affinity-chromatography using	
	epitope-specific immunogen.	
Clonality	Polyclonal	
Concentration	1 mg/ml	
Observed band	190kD	
Human Gene ID	2909	
Human Swiss-Prot Number	Q9NRY4	
Alternative Names	ARHGAP35; GRF1; GRLF1; KIAA1722; Rho	
	GTPase-activating protein 35; Glucocorticoid	



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



Background

293T

1/0---70--55---40--35---25---

15----

1.2

1

0.8

0.4

0.2

0D 450 nm

p-GRF-1 (Y1087)

receptor DNA-binding factor 1; Glucocorticoid receptor repression factor 1; GRF-1; Rho GAP p190A; p190-A

The human glucocorticoid receptor DNA binding factor, which associates with the promoter region of the glucocorticoid receptor gene (hGR gene), is a repressor of glucocorticoid receptor transcription. The amino acid sequence deduced from the cDNA sequences show the presence of three sequence motifs characteristic of a zinc finger and one motif suggestive of a leucine zipper in which 1 cysteine is found instead of all leucines. The GRLF1 enhances the homologous down-regulation of wild-type hGR gene expression. Biochemical analysis suggests that GRLF1 interaction is sequence specific and that transcriptional efficacy of GRLF1 is regulated through its interaction with specific sequence motif. The level of expression is regulated by glucocorticoids. [provided by RefSeq, Jul 2008],

Western Blot analysis of 293T cells using Phospho-GRF-1 (Y1087) Polyclonal Antibody diluted at 1:1000

Enzyme-Linked Immunosorbent Assay (Phospho-ELISA) for Immunogen Phosphopeptide (Phospho-left) and Non-Phosphopeptide (Phospho-right), using GRF-1 (Phospho-Tyr1087) Antibody

+86-27-59760950

Non-Phospho

Phosphopeptide

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C

ELK Biotechnology		
	GRF-1	Western blot analysis of lysates from K562 cells, using phospho-GRF-1 (Phospho-Tyr1087) antibody.
	-95	
	-72 -55	



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road,Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C