

## mGluR-8 rabbit pAb

Cat No.: ES5680

For research use only

## Overview

Product Name mGluR-8 rabbit pAb

Host species Rabbit

Applications WB;IHC;IF;ELISA Species Cross-Reactivity Human;Mouse;Rat

**Recommended dilutions** Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications. The antiserum was produced against synthesized

Immunogen The antiserum was produced against synthesized

peptide derived from human mGluR8. AA

range:841-890

**Specificity** mGluR-8 Polyclonal Antibody detects endogenous

levels of mGluR-8 protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

**Store at -20°C.** Avoid repeated freeze-thaw cycles.

Protein Name Metabotropic glutamate receptor 8

Gene Name GRM8

Cellular localizationCell membrane; Multi-pass membrane protein.PurificationThe antibody was affinity-purified from rabbit antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 100kD
Human Gene ID 2918
Human Swiss-Prot Number 000222

Alternative Names GRM8; GPRC1H; MGLUR8; Metabotropic glutamate

receptor 8; mGluR8

**Background** glutamate metabotropic receptor 8(GRM8) Homo

sapiens L-glutamate is the major excitatory neurotransmitter in the central nervous system and

activates both ionotropic and metabotropic

glutamate receptors. Glutamatergic




+86-27-59760950 ELKbio@ELKbiotech.com


www.elkbiotech.com




neurotransmission is involved in most aspects of normal brain function and can be perturbed in many neuropathologic conditions. The metabotropic glutamate receptors are a family of G protein-coupled receptors, that have been divided into 3 groups on the basis of sequence homology, putative signal transduction mechanisms, and pharmacologic properties. Group I includes GRM1 and GRM5 and these receptors have been shown to activate phospholipase C. Group II includes GRM2 and GRM3 while Group III includes GRM4, GRM6, GRM7 and GRM8. Group II and III receptors are linked to the inhibition of the cyclic AMP cascade but differ in their agonist selectivities. Alternatively spliced transcript variants encoding different isoforms have been d

Immunohistochemistry analysis of paraffin-embedded human brain tissue, using mGluR8 Antibody. The picture on the right is blocked with the synthesized peptide.







+86-27-59760950

Western blot analysis of lysates from mouse brain, using mGluR8 Antibody. The lane on the right is blocked with the synthesized peptide.

