

HXK III rabbit pAb

Cat No.:ES5740

For research use only

Overview

Product Name	HXK III rabbit pAb
Host species	Rabbit
Applications	WB;IHC;IF;ELISA
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	Western Blot: 1/500 - 1/2000.
	Immunohistochemistry: 1/100 - 1/300. ELISA:
	1/5000. Not yet tested in other applications.
Immunogen	The antiserum was produced against synthesized
	peptide derived from human Hexokinase-3. AA
	range:811-860
Specificity	HXK III Polyclonal Antibody detects endogenous
	levels of HXK III protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Hexokinase-3
Gene Name	HK3
Cellular localization	cell,mitochondrion,cytosol,
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	99kD
Human Gene ID	3101
Human Swiss-Prot Number	P52790
Alternative Names	HK3; Hexokinase-3; Hexokinase type III; HK III
Background	Hexokinases phosphorylate glucose to produce
	glucose-6-phosphate, the first step in most glucose
	metabolism pathways. This gene encodes
	hexokinase 3. Similar to hexokinases 1 and 2, this
	allosteric enzyme is inhibited by its product
	glucose-6-phosphate. [provided by RefSeq, Apr





+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C





Immunohistochemical analysis of paraffin-embedded human spleen tissue. 1,primary Antibody was diluted at 1:200(4° overnight). 2, Sodium citrate pH 6.0 was used for antigen retrieval(>98°C,20min). 3,Secondary antibody was diluted at 1:200



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C