

## HMG-17 rabbit pAb

Cat No.: ES5746

For research use only

## Overview

Product Name HMG-17 rabbit pAb

Host species Rabbit
Applications IHC;IF;ELISA

Species Cross-Reactivity Human; Mouse; Rat

**Recommended dilutions** Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.

Immunogen The antiserum was produced against synthesized

peptide derived from human HMG17. AA range:1-50

Specificity HMG-17 Polyclonal Antibody detects endogenous

levels of HMG-17 protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Non-histone chromosomal protein HMG-17

Gene Name HMGN2

**Cellular localization** Nucleus . Cytoplasm . Cytoplasmic enrichment upon

phosphorylation.

**Purification** The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal Concentration 1 mg/ml

**Observed band** 

Human Gene ID 3151 Human Swiss-Prot Number P05204

Alternative Names HMGN2; HMG17; Non-histone chromosomal protein

HMG-17; High mobility group nucleosome-binding

domain-containing protein 2

Background high mobility group nucleosomal binding domain

2(HMGN2) Homo sapiens The protein encoded

by this gene binds nucleosomal DNA and is

associated with transcriptionally active chromatin.




+86-27-59760950 ELKbio@ELKbiotech.com www.elkbiotech.com




Along with a similar protein, HMGN1, the encoded protein may help maintain an open chromatin configuration around transcribable genes. The protein has also been found to have antimicrobial activity against bacteria, viruses and fungi. [provided by RefSeq, Oct 2014],

Immunofluorescence analysis of HeLa cells, using HMG17 Antibody. The picture on the right is blocked with the synthesized peptide.



Immunohistochemistry analysis of paraffin-embedded human brain tissue, using HMG17 Antibody. The picture on the right is blocked with the synthesized peptide.



+86-27-59760950

