

ILK (phospho Ser246) rabbit pAb

Cat No.: ES5902

For research use only

Overview

Product Name ILK (phospho Ser246) rabbit pAb

Host species Rabbit
Applications IHC;IF;ELISA

Species Cross-Reactivity Human; Mouse; Rat

Recommended dilutions Immunohistochemistry: 1/100 - 1/300. ELISA:

1/5000. Not yet tested in other applications.

Immunogen The antiserum was produced against synthesized

peptide derived from human ILK around the

phosphorylation site of Ser246. AA range:212-261

Specificity Phospho-ILK (S246) Polyclonal Antibody detects

endogenous levels of ILK protein only when

phosphorylated at S246.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Integrin-linked protein kinase

Gene Name ILK

Cellular localization Cell junction, focal adhesion . Cell membrane;

Peripheral membrane protein; Cytoplasmic side .

Cell projection, lamellipodium. Cytoplasm,

myofibril, sarcomere.

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal Concentration 1 mg/ml

Observed band

Human Gene ID 3611 Human Swiss-Prot Number Q13418

Alternative Names ILK; ILK1; ILK2; Integrin-linked protein kinase; 59 kDa

serine/threonine-protein kinase; ILK-1; ILK-2; p59ILK

Background This gene encodes a protein with a kinase-like

domain and four ankyrin-like repeats. The encoded




+86-27-59760950 ELKbio@ELKbiotech.com www.elkbiotech.com



protein associates at the cell membrane with the cytoplasmic domain of beta integrins, where it regulates integrin-mediated signal transduction. Activity of this protein is important in the epithelial to mesenchymal transition, and over-expression of this gene is implicated in tumor growth and metastasis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jun 2013],

Immunohistochemistry analysis of paraffin-embedded human brain, using ILK (Phospho-Ser246) Antibody. The picture on the right is blocked with the phospho peptide.



+86-27-59760950

