

FAM111B rabbit pAb

Cat No.: ES5990

For research use only

Overview

Product Name FAM111B rabbit pAb

Host species Rabbit
Applications IHC;IF;ELISA

Species Cross-Reactivity Human; Rat; Mouse;

Recommended dilutions Immunohistochemistry: 1/100 - 1/200. ELISA:

1/5000. Not yet tested in other applications.

Immunogen The antiserum was produced against synthesized

peptide derived from human F111B. AA

range:281-330

Specificity FAM111B Polyclonal Antibody detects endogenous

levels of FAM111B protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Protein FAM111B

Gene Name FAM111B

Cellular localization

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal Concentration 1 mg/ml

Observed band

Human Gene ID 374393 Human Swiss-Prot Number Q6SJ93

Alternative Names FAM111B; CANP; Protein FAM111B;

Cancer-associated nucleoprotein

Background This gene encodes a protein with a trypsin-like

cysteine/serine peptidase domain in the C-terminus.

Mutations in this gene are associated with an autosomal dominant form of hereditary fibrosing poikiloderma (HFP). Affected individuals display mottled pigmentation, telangiectasia, epidermal



+86-27-59760950 ELKbio@ELKbiotech.com ww




atrophy, tendon contractures, and progressive pulmonary fibrosis. Alternative splicing results in multiple transcript variants encoding distinct isoforms. A paralog of this gene which also has a trypsin-like peptidase domain, FAM111A, is located only 16 kb from this gene on human chromosome 11q12.1. [provided by RefSeq, Apr 2014],

Immunohistochemical analysis of paraffin-embedded

Human brain. Antibody was diluted at 1:100(4° overnight). High-pressure and temperature

Tris-EDTA,pH8.0 was used for antigen retrieval. Negetive contrl (right) obtaned from antibody was pre-absorbed by i



+86-27-59760950

