

## KV4.1 rabbit pAb

## Cat No.:ES5996

For research use only

## Overview

Product Name	KV4.1 rabbit pAb
Host species	Rabbit
Applications	WB;IHC;IF;ELISA
Species Cross-Reactivity	Human; Mouse
<b>Recommended dilutions</b>	Western Blot: 1/500 - 1/2000.
	Immunohistochemistry: 1/100 - 1/300. ELISA:
	1/20000. Not yet tested in other applications.
Immunogen	The antiserum was produced against synthesized
	peptide derived from human KCND1. AA
	range:558-607
Specificity	KV4.1 Polyclonal Antibody detects endogenous
	levels of KV4.1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Potassium voltage-gated channel subfamily D
	member 1
Gene Name	KCND1
Cellular localization	Membrane; Multi-pass membrane protein. Cell
	projection, dendrite .
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	70kD
Human Gene ID	3750
Human Swiss-Prot Number	Q9NSA2
Alternative Names	KCND1; Potassium voltage-gated channel subfamily
	D member 1; Voltage-gated potassium channel
	subunit Kv4.1
Background	This gene encodes a multipass membrane protein
	that comprises the pore subunit of the



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com


23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



voltage-gated A-type potassium channel, which functions in the repolarization of membrane action potentials. Activity of voltage-gated potassium channels is important in a number of physiological processes, among them the regulation of neurotransmitter release, heart rate, insulin secretion, and smooth muscle contraction. [provided by RefSeq, Aug 2013],

Western Blot analysis of 3T3 cells using KV4.1 Polyclonal Antibody

Western blot analysis of lysates from NIH/3T3 cells, using KCND1 Antibody. The lane on the right is blocked with the synthesized peptide.





www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C