

## MDMX rabbit pAb

Cat No.:ES6214

For research use only

## Overview

Product Name MDMX rabbit pAb

Host species Rabbit

Applications WB;IHC;IF;ELISA Species Cross-Reactivity Human;Mouse;Rat

**Recommended dilutions** Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.

Immunogen The antiserum was produced against synthesized

peptide derived from human MDM4. AA

range:336-385

**Specificity** MDMX Polyclonal Antibody detects endogenous

levels of MDMX protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

**Storage** Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Protein Mdm4

**Gene Name** MDM4 **Cellular localization** Nucleus.

**Purification** The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 45kD
Human Gene ID 4194
Human Swiss-Prot Number O15151

Alternative Names MDM4; MDMX; Protein Mdm4; Double minute 4

protein; Mdm2-like p53-binding protein; Protein

Mdmx; p53-binding protein Mdm4

**Background** This gene encodes a nuclear protein that contains a

p53 binding domain at the N-terminus and a RING

finger domain at the C-terminus, and shows

structural similarity to p53-binding protein MDM2.




+86-27-59760950 ELKbio@ELKbiotech.com www.elkbiotech.co



Both proteins bind the p53 tumor suppressor protein and inhibit its activity, and have been shown to be overexpressed in a variety of human cancers. However, unlike MDM2 which degrades p53, this protein inhibits p53 by binding its transcriptional activation domain. This protein also interacts with MDM2 protein via the RING finger domain, and inhibits the latter's degradation. So this protein can reverse MDM2-targeted degradation of p53, while maintaining suppression of p53 transactivation and apoptotic functions. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene. [provided by RefSeq, Feb 2011],




Western Blot analysis of 293T cells using MDMX Polyclonal Antibody diluted at 1:1000 cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003,Inventbiotech,MN,USA).



Immunohistochemistry analysis of paraffin-embedded human colon carcinoma tissue, using MDM4 Antibody. The picture on the right is blocked with the synthesized peptide.







Western blot analysis of lysates from rat muscle cells, using MDM4 Antibody. The lane on the right is blocked with the synthesized peptide.

