

## ANKRD20 rabbit pAb

Cat No.: ES6272

For research use only

## Overview

Product Name ANKRD20 rabbit pAb

Host species Rabbit
Applications IHC;IF;ELISA

Species Cross-Reactivity Human; Rat; Mouse;

**Recommended dilutions** Immunohistochemistry: 1/100 - 1/300. ELISA:

1/10000. Not yet tested in other applications.

Immunogen The antiserum was produced against synthesized

peptide derived from human ANKRD20A3. AA

range:551-600

Specificity ANKRD20 Polyclonal Antibody detects endogenous

levels of ANKRD20 protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Ankyrin repeat domain-containing protein 20A3

Gene Name ANKRD20A3

**Cellular localization** 

**Purification** The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 95kD
Human Gene ID 441425
Human Swiss-Prot Number Q5VUR7

Alternative Names ANKRD20A3; Ankyrin repeat domain-containing

protein 20A3


**Background** similarity:Contains 5 ANK repeats.,tissue

specificity: Widely expressed; not found in kidney,

nor in liver.,







Immunohistochemistry analysis of paraffin-embedded human brain, using ANKRD20A3 Antibody. The picture on the right is blocked with the synthesized peptide.



+86-27-59760950