

## RFPL4B rabbit pAb

Cat No.:ES6280

For research use only

## Overview

Product Name RFPL4B rabbit pAb

Host species Rabbit
Applications IHC;IF;ELISA

**Species Cross-Reactivity** Human;Rat;Mouse;

**Recommended dilutions** Immunohistochemistry: 1/100 - 1/300. ELISA:

1/40000. Not yet tested in other applications.

Immunogen The antiserum was produced against synthesized

peptide derived from human RFPL4B. AA

range:100-149

Specificity RFPL4B Polyclonal Antibody detects endogenous

levels of RFPL4B protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

**Storage** Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Ret finger protein-like 4B

Gene Name RFPL4B

**Cellular localization** 

**Purification** The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal Concentration 1 mg/ml

**Observed band** 

Human Gene ID 442247 Human Swiss-Prot Number Q6ZWI9

Alternative Names RFPL4B; RNF211; Ret finger protein-like 4B; RING

finger protein 211

Background similarity: Contains 1 B30.2/SPRY

domain., similarity: Contains 1 RING-type zinc finger.,



+86-27-59760950 ELKbio@ELKbiotech.com

www.elkbiotech.com





Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 30min).



+86-27-59760950