

DHHC-9 rabbit pAb

Cat No.: ES6497

For research use only

Overview

Product Name DHHC-9 rabbit pAb

Host species Rabbit

Applications WB;IHC;IF;ELISA Species Cross-Reactivity Human;Mouse;Rat

Recommended dilutions Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.

Immunogen The antiserum was produced against synthesized

peptide derived from human ZDHHC9. AA

range:315-364

Specificity DHHC-9 Polyclonal Antibody detects endogenous

levels of DHHC-9 protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Palmitoyltransferase ZDHHC9

Gene Name ZDHHC9

Cellular localization Endoplasmic reticulum membrane; Multi-pass

membrane protein . Golgi apparatus membrane ;

Multi-pass membrane protein .

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 41kD
Human Gene ID 51114
Human Swiss-Prot Number Q9Y397

Alternative Names ZDHHC9; CXorf11; ZDHHC10; ZNF379; ZNF380;


CGI-89; Palmitoyltransferase ZDHHC9; Zinc finger DHHC domain-containing protein 9; DHHC-9; DHHC9; Zinc finger protein 379; Zinc finger protein

380






Background



This gene encodes an integral membrane protein that is a member of the zinc finger DHHC domain-containing protein family. The encoded protein forms a complex with golgin subfamily A member 7 and functions as a palmitoyltransferase. This protein specifically palmitoylates HRAS and NRAS. Mutations in this gene are associated with X-linked mental retardation. Alternate splicing results in multiple transcript variants that encode the same protein.[provided by RefSeq, May 2010],

Western Blot analysis of 293T cells using DHHC-9 Polyclonal Antibody diluted at 1:1000




+86-27-59760950

Immunohistochemistry analysis of paraffin-embedded human brain tissue, using ZDHHC9 Antibody. The picture on the right is blocked with the synthesized peptide.



23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei, P.R.C





Western blot analysis of lysates from COLO cells, using ZDHHC9 Antibody. The lane on the right is blocked with the synthesized peptide.

