

PEG3 rabbit pAb

Cat No.:ES6567

For research use only

Overview

Product Name	PEG3 rabbit pAb
Host species	Rabbit
Applications	WB;IHC
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:500-2000;IHC-p 1:50-300
Immunogen	The antiserum was produced against synthesized
	peptide derived from human PEG3. AA
	range:1031-1080
Specificity	PEG3 Polyclonal Antibody detects endogenous levels
	of PEG3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Paternally-expressed gene 3 protein
Gene Name	PEG3
Cellular localization	Nucleus . Cytoplasm .
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	181kD
Human Gene ID	5178
Human Swiss-Prot Number	Q9GZU2
Alternative Names	PEG3; KIAA0287; ZSCAN24; Paternally-expressed
	gene 3 protein; Zinc finger and SCAN
	domain-containing protein 24
Background	In human, ZIM2 and PEG3 are treated as two
	distinct genes though they share multiple 5' exons
	and a common promoter and both genes are
	paternally expressed (PMID:15203203). Alternative
	splicing events connect their shared 5' exons either
_	with the remaining 4 exons unique to ZIM2, or with



+86-27-59760950

ELKbio@ELKbiotech.com


www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C




the remaining 2 exons unique to PEG3. In contrast, in other mammals ZIM2 does not undergo imprinting and, in mouse, cow, and likely other mammals as well, the ZIM2 and PEG3 genes do not share exons. Human PEG3 protein belongs to the Kruppel C2H2-type zinc finger protein family. PEG3 may play a role in cell proliferation and p53-mediated apoptosis. PEG3 has also shown tumor suppressor activity and tumorigenesis in glioma and ovarian cells. Alternative splicing of this PEG3 gene results in multiple transcript variants encoding distinct isoforms. [provided by RefSeq, Sep 2009],

Western Blot analysis of A549 cells using PEG3 Polyclonal Antibody



Immunohistochemical analysis of paraffin-embedded human cervical carcinoma. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).





+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C