

## Girdin (phospho Ser1417) rabbit pAb

Cat No.: ES6748

For research use only

## Overview

Product Name Girdin (phospho Ser1417) rabbit pAb

Host speciesRabbitApplicationsIHC;IF;ELISASpecies Cross-ReactivityHuman;Mouse

**Recommended dilutions** Immunohistochemistry: 1/100 - 1/300. ELISA:

1/5000. Not yet tested in other applications.

Immunogen The antiserum was produced against synthesized

peptide derived from human Girdin around the

phosphorylation site of Ser1417. AA

range:1383-1432

**Specificity** Phospho-Girdin (S1417) Polyclonal Antibody detects

endogenous levels of Girdin protein only when

phosphorylated at S1417.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

**Storage** Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Girdin
Gene Name CCDC88A

Cellular localization Cell membrane; Peripheral membrane protein.

Cytoplasm, cytosol . Cytoplasmic vesicle . Cell

projection, lamellipodium. Cytoplasm, cytoskeleton,

cilium basal body . Cytoplasm, cytoskeleton, microtubule organizing center, centrosome,

centriole. Localize

**Purification** The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal Concentration 1 mg/ml

**Observed band** 

Human Gene ID 55704 Human Swiss-Prot Number Q3V6T2

Alternative Names CCDC88A; APE; GRDN; KIAA1212; Girdin; Akt




+86-27-59760950 ELKbio@ELKbiotech.com

www.elkbiotech.com




**Background** 

phosphorylation enhancer; APE; Coiled-coil domain-containing protein 88A; G alpha-interacting vesicle-associated protein; GIV; Girders of actin filament; Hook-related protein 1; HkRP1 This gene encodes a member of the Girdin family of coiled-coil domain containing proteins. The encoded protein is an actin-binding protein that is activated by the serine/threonine kinase Akt and plays a role in cytoskeleton remodeling and cell migration. The encoded protein also enhances Akt signaling by mediating phosphoinositide 3-kinase (PI3K)-dependent activation of Akt by growth factor receptor tyrosine kinases and G protein-coupled receptors. Increased expression of this gene and phosphorylation of the encoded protein may play a role in cancer metastasis. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Dec 2011],



Enzyme-Linked Immunosorbent Assay (Phospho-ELISA) for Immunogen Phosphopeptide (Phospho-left) and Non-Phosphopeptide (Phospho-right), using Girdin (Phospho-Ser1417) Antibody



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma, using Girdin (Phospho-Ser1417) Antibody. The picture on the right is blocked with the phospho peptide.

