

PKC δ (phospho Thr507) rabbit pAb

Cat No.:ES6772

For research use only

Overview

Product Name PKC δ (phospho Thr507) rabbit pAb

Host species Rabbit

Applications WB;IHC;IF;ELISA Species Cross-Reactivity Human;Mouse;Rat

Recommended dilutions Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.

Immunogen The antiserum was produced against synthesized

peptide derived from human PKC delta around the phosphorylation site of Thr505. AA range:471-520

Specificity Phospho-PKC δ (T507) Polyclonal Antibody detects

endogenous levels of PKC $\boldsymbol{\delta}$ protein only when

phosphorylated at T507.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Protein kinase C delta type

Gene Name PRKCD

Cellular localization Cytoplasm . Cytoplasm, perinuclear region .

Nucleus . Cell membrane ; Peripheral membrane protein . Mitochondrion . Endomembrane system . Translocates to the mitochondria upon apoptotic stimulation. Upon activation, translocates to the plasma membrane followed by partial location to

the endolysosomes (PubMed:17303575). .

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 78kD
Human Gene ID 5580
Human Swiss-Prot Number Q05655



+86-27-59760950 ELKbio@ELKbiotech.com www.elkbiotech.com




Alternative Names

Background


PRKCD; Protein kinase C delta type; Tyrosine-protein kinase PRKCD; nPKC-delta

Protein kinase C (PKC) is a family of serine- and threonine-specific protein kinases that can be activated by calcium and the second messenger diacylglycerol. PKC family members phosphorylate a wide variety of protein targets and are known to be involved in diverse cellular signaling pathways. PKC family members also serve as major receptors for phorbol esters, a class of tumor promoters. Each member of the PKC family has a specific expression profile and is believed to play distinct roles in cells. The protein encoded by this gene is one of the PKC family members. Studies both in human and mice demonstrate that this kinase is involved in B cell signaling and in the regulation of growth, apoptosis, and differentiation of a variety of cell types. Alternatively spliced transcript variants encoding the same protein have been observed. [provided by

RefSeq, Jul 2008],



Western Blot analysis of 3T3 cells using Phospho-PKC δ (T507) Polyclonal Antibody



Enzyme-Linked Immunosorbent Assay (Phospho-ELISA) for Immunogen Phosphopeptide (Phospho-left) and Non-Phosphopeptide (Phospho-right), using PKC delta (Phospho-Thr505) Antibody




ELKbio@ELKbiotech.com


www.elkbiotech.com







Immunohistochemistry analysis of paraffin-embedded human breast carcinoma, using PKC delta (Phospho-Thr505) Antibody. The picture on the right is blocked with the phospho peptide.



Western blot analysis of lysates from NIH/3T3 cells treated with UV 15', using PKC delta (Phospho-Thr505) Antibody. The lane on the right is blocked with the phospho peptide.

