

PKC δ (phospho Tyr64) rabbit pAb

Cat No.: ES6773

For research use only

Overview

Product Name PKC δ (phospho Tyr64) rabbit pAb

Host species Rabbit WB;IHC **Applications**

Species Cross-Reactivity Human; Mouse; Rat

Recommended dilutions WB 1:500-2000;IHC-p 1:50-300

Immunogen The antiserum was produced against synthesized

> peptide derived from human PKCD around the phosphorylation site of Tyr64. AA range:30-79 Phospho-PKC δ (Y64) Polyclonal Antibody detects

Specificity endogenous levels of PKC δ protein only when

phosphorylated at Y64.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Store at -20°C. Avoid repeated freeze-thaw cycles. Storage

Protein Name Protein kinase C delta type

Gene Name PRKCD

Cellular localization Cytoplasm . Cytoplasm, perinuclear region .

> Nucleus . Cell membrane ; Peripheral membrane protein. Mitochondrion. Endomembrane system. Translocates to the mitochondria upon apoptotic stimulation. Upon activation, translocates to the

> > www.elkbiotech.com

plasma membrane fol

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal Concentration 1 mg/ml **Observed band** 77kD **Human Gene ID** 5580 **Human Swiss-Prot Number** Q05655

Alternative Names PRKCD; Protein kinase C delta type; Tyrosine-protein

kinase PRKCD; nPKC-delta

Background Protein kinase C (PKC) is a family of serine- and



+86-27-59760950 ELKbio@ELKbiotech.com




threonine-specific protein kinases that can be activated by calcium and the second messenger diacylglycerol. PKC family members phosphorylate a wide variety of protein targets and are known to be involved in diverse cellular signaling pathways. PKC family members also serve as major receptors for phorbol esters, a class of tumor promoters. Each member of the PKC family has a specific expression profile and is believed to play distinct roles in cells. The protein encoded by this gene is one of the PKC family members. Studies both in human and mice demonstrate that this kinase is involved in B cell signaling and in the regulation of growth, apoptosis, and differentiation of a variety of cell types. Alternatively spliced transcript variants encoding the same protein have been observed. [provided by RefSeq, Jul 2008],

0.500
0.450
0.400
0.350
0.300
0.250
0.200
0.150
0.100
0.050
0.000

phosphopeptide non-phosphopeptide

Enzyme-Linked Immunosorbent Assay (Phospho-ELISA) for Immunogen Phosphopeptide (Phospho-left) and Non-Phosphopeptide (Phospho-right), using PKCD (Phospho-Tyr64) Antibody



+86-27-59760950


Western blot analysis of lysates from HepG2 cells treated with PMA 125ng/ml 30', using PKCD (Phospho-Tyr64) Antibody. The lane on the right is blocked with the phospho peptide.



ELKbio@ELKbiotech.com

www.elkbiotech.com





Immunohistochemical analysis of paraffin-embedded human uterus. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).



+86-27-59760950