

Presenilin 2 (phospho Ser330) rabbit pAb

Cat No.: ES6849

For research use only

Overview

Product Name Presenilin 2 (phospho Ser330) rabbit pAb

Host species Rabbit

Applications WB;IHC;IF;ELISA Species Cross-Reactivity Human;Mouse;Rat

Recommended dilutions Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.

Immunogen The antiserum was produced against synthesized

peptide derived from human Presenilin 2 around the phosphorylation site of Ser330. AA range:296-345

Specificity Phospho-Presenilin 2 (S330) Polyclonal Antibody

detects endogenous levels of Presenilin 2 protein

only when phosphorylated at S330.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Presenilin-2
Gene Name PSEN2

Cellular localization Endoplasmic reticulum membrane ; Multi-pass

membrane protein . Golgi apparatus membrane ;

Multi-pass membrane protein.

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 50kD
Human Gene ID 5664
Human Swiss-Prot Number P49810

Alternative Names PSEN2; AD4; PS2; PSNL2; STM2; Presenilin-2; PS-2;

AD3LP; AD5; E5-1; STM-2

Background Alzheimer's disease (AD) patients with an inherited

form of the disease carry mutations in the presenilin




+86-27-59760950 ELKbio@ELKbiotech.com www.elkbiotech.co




proteins (PSEN1 or PSEN2) or the amyloid precursor protein (APP). These disease-linked mutations result in increased production of the longer form of amyloid-beta (main component of amyloid deposits found in AD brains). Presenilins are postulated to regulate APP processing through their effects on gamma-secretase, an enzyme that cleaves APP. Also, it is thought that the presenilins are involved in the cleavage of the Notch receptor such that, they either directly regulate gamma-secretase activity, or themselves act are protease enzymes. Two alternatively spliced transcript variants encoding different isoforms of PSEN2 have been identified. [provided by RefSeq, Jul 2008],

Immunohistochemical analysis of paraffin-embedded Human breast cancer. Antibody was diluted at 1:100(4° overnight). High-pressure and temperature Tris-EDTA,pH8.0 was used for antigen retrieval. Negetive contrl (right) obtaned from antibody was pre-absor




Immunohistochemical analysis of paraffin-embedded Human colon cancer. Antibody was diluted at 1:100(4° overnight). High-pressure and temperature Tris-EDTA,pH8.0 was used for antigen retrieval. Negetive contrl (right) obtaned from antibody was pre-absorb



+86-27-59760950







Immunohistochemical analysis of paraffin-embedded Human colon cancer. Antibody was diluted at 1:100(4° overnight). High-pressure and temperature Tris-EDTA,pH8.0 was used for antigen retrieval. Negetive contrl (right) obtaned from antibody was pre-absorb



+86-27-59760950