

## CysLTR2 rabbit pAb

Cat No.: ES6868

For research use only

## Overview

Product Name CysLTR2 rabbit pAb

Host species Rabbit
Applications WB;IF;ELISA

Species Cross-Reactivity Human; Rat; Mouse;

**Recommended dilutions** Western Blot: 1/500 - 1/2000. Immunofluorescence:

1/200 - 1/1000. ELISA: 1/20000. Not yet tested in

other applications.

Immunogen The antiserum was produced against synthesized

peptide derived from human CYSLTR2. AA

range:281-330

**Specificity** CysLTR2 Polyclonal Antibody detects endogenous

levels of CysLTR2 protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

**Storage** Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Cysteinyl leukotriene receptor 2

Gene Name CYSLTR2

Cellular localizationCell membrane; Multi-pass membrane protein.PurificationThe antibody was affinity-purified from rabbit antiserum by affinity-chromatography using

epitope-specific immunogen.

ClonalityPolyclonalConcentration1 mg/mlObserved band40kDHuman Gene ID57105Human Swiss-Prot NumberQ9NS75

Alternative Names CYSLTR2; CYSLT2; CYSLT2R; PSEC0146; Cysteinyl

leukotriene receptor 2; CysLTR2; G-protein coupled receptor GPCR21; hGPCR21; G-protein coupled

receptor HG57; HPN321

Background The cysteinyl leukotrienes LTC4, LTD4, and LTE4 are


important mediators of human bronchial asthma. Pharmacologic studies have determined that






cysteinyl leukotrienes activate at least 2 receptors, the protein encoded by this gene and CYSLTR1. This encoded receptor is a member of the superfamily of G protein-coupled receptors. It seems to play a major role in endocrine and cardiovascular systems. [provided by RefSeq, Jul 2008],

Western Blot analysis of COLO cells using CysLTR2 Polyclonal Antibody



Western blot analysis of lysates from COLO cells, using CYSLTR2 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HT-29 cells using CYSLTR2 antibody.

